

INSTRUCTION MANUAL

STEFAN'S CONSTANT KIT



MITTAL ENTERPRISES

2151/T-7C, New Patel Nagar, New Delhi-110008

Ph. : 25702784 (Off.) 25120261 (Res.); Fax : 25120261

Mobile: 9810681132

E-mail : mittalenterprises@bol.net.in

STEFAN'S CONSTANT KIT

INTRODUCTION

According to Stefan-Boltzmann's law the heat energy E_1 radiated per unit area per second by a body at T° Kelvin surrounded by another body at lower temperature T_0° Kelvin is directly proportional to the difference in fourth power of the absolute temperature i.e.

$$E_1 \propto (T^4 - T_0^4)$$

$$\text{or} \quad E_1 = \sigma \cdot (T^4 - T_0^4)$$

where σ is called Stefan's Constant. Its standard value is $5.67 \times 10^{-8} \text{ W m}^{-2} \text{ K}^{-4}$. If the area of surface is A , then the heat energy radiated from the surface becomes

$$E_1 = \sigma \cdot A \cdot (T^4 - T_0^4) \quad \dots(1)$$

THEORY

In this kit, the electric energy is used to heat up the two copper circular discs to a particular steady temperature $T^\circ \text{ K}$. The steady temperature is measured by thermometers placed in contact with the copper discs. As the steady temperature of the discs is larger than surrounding temperature $T_0^\circ \text{ K}$ of the atmosphere, the disc will radiate heat energy which can be calculated by Stefan's Boltzmann's Law or eqn. (1).

As the temperature of the discs reaches steady state, heat energy, E_1 , radiated by discs per second will be equal to the electrical energy, E_2 , consumed by heating element per second and is given by

$$E_2 = V \cdot I \quad \dots(2)$$

where V = Potential difference across the heating element and
 I = Current flowing in the heating element

From equations (1) & (2)

$$V \cdot I = \sigma \cdot A \cdot (T^4 - T_0^4)$$

$$\text{or} \quad \sigma = \frac{V \cdot I}{A \cdot (T^4 - T_0^4)} \quad \dots(3)$$

V and I can be measured with the help of Voltmeter and Ammeter installed in the kit. The steady temperature (T °K) of copper discs can be read from the thermometers and the area A of the disc can be calculated by measuring the radius of the disc.

By substituting the values of V, I, A, T and T_0 in equation (3) one can calculate the Stefan's constant.

The features of the kit are:

- i) There is practically no loss of heat energy from heating element by conduction, convection or radiation. Therefore, entire heat energy of the heating element is transferred to discs.
- ii) As the efficiency of the heating element is nearly equal to one, almost entire electrical energy supplied to the heater will be converted into heat energy.

PROCEDURE

- i) Note down the surrounding room temperature (T °K) with the help of the thermometers.
- ii) Put ON the switch on the kit.
- iii) Place the discs horizontally on the table.
- iv) Connect the heating element of discs to the Kit terminals marked "H".
- v) Adjust with the help of potentiometer, the potential difference (V) and current (I) for heating element for proper reading in the meters.
- vi) Place the thermometers vertically in the holes provided in the discs.
- vii) Observe the readings of the thermometers till the readings become steady. Take the average of the readings.

- viii) The several sets of observations are taken by varying the potential difference (V) & current (I) with the help of potentiometer.
- ix) The area of the plate can be calculated by measuring the radius of the disc.

SAMPLE OBSERVATIONS

The area of the copper plates = 306.8056 sq. cm.

Least Count of Thermometer = 1 °C.

S.No.	Voltage	Current	Room Temp. T °K	Mean Temp. of Thermometers T °K
1	20	0.25	300	328
2	30	0.35	300	395
3	40	0.47	300	376
4	50	0.58	300	404

SAMPLE CALCULATIONS

$$\sigma = \frac{V \cdot I}{A \cdot (T^4 - T_0^4)} \text{ W cm}^{-2} \text{ K}^{-4}$$

$$\sigma = \frac{V \cdot I \times 10^7}{A \cdot (T^4 - T_0^4)} \text{ ergs cm}^{-2} \text{ sec}^{-1} \text{ K}^{-4}$$

Thus substituting values from table :

$$\begin{aligned} 1) \quad \sigma &= \frac{20 \times 0.25 \times 10^7}{306.8056 \cdot (328^4 - 300^4)} \\ &= \frac{20 \times 0.25 \times 10^7}{306.8056 \times 3.482 \times 10^9} \\ &= 4.679 \times 10^{-5} \text{ erg cm}^{-2} \text{ sec}^{-1} \text{ K}^{-4} \end{aligned}$$

$$\begin{aligned} 2) \quad \sigma &= \frac{30 \times 0.35 \times 10^7}{306.8056 \cdot (345^4 - 300^4)} \\ &= \frac{20 \times 0.35 \times 10^7}{306.8056 \times 6.072 \times 10^9} \\ &= 5.637 \times 10^{-5} \text{ erg cm}^{-2} \text{ sec}^{-1} \text{ K}^{-4} \end{aligned}$$

$$\begin{aligned} 3) \quad \sigma &= \frac{40 \times 0.47 \times 10^7}{306.8056 \cdot (376^4 - 300^4)} \\ &= \frac{40 \times 0.47 \times 10^7}{306.8056 \times 11.892 \times 10^9} \\ &= 5.141 \times 10^{-5} \text{ erg cm}^{-2} \text{ sec}^{-1} \text{ K}^{-4} \end{aligned}$$

$$\begin{aligned} 4) \quad \sigma &= \frac{50 \times 0.58 \times 10^7}{306.8056 \cdot (404^4 - 300^4)} \\ &= \frac{50 \times 0.58 \times 10^7}{306.8056 \times 18.102 \times 10^9} \\ &= 5.223 \times 10^{-5} \text{ erg cm}^{-2} \text{ sec}^{-1} \text{ K}^{-4} \end{aligned}$$

$$\text{Average} \quad \sigma = 5.17 \times 10^{-5} \text{ erg cm}^{-2} \text{ sec}^{-1} \text{ K}^{-4}$$

Note: All above measurements and readings varies from instrument to instrument due to change in dimensions.

- viii) The several sets of observations are taken by varying the potential difference (V) & current (I) with the help of potentiometer.
- ix) The area of the plate can be calculated by measuring the radius of the disc.

SAMPLE OBSERVATIONS

The area of the copper plates = 306.8056 sq. cm.

Least Count of Thermometer = 1 °C.

S.No.	Voltage	Current	Room Temp. T °K	Mean Temp. of Thermometers T °K
1	20	0.25	300	328
2	30	0.35	300	395
3	40	0.47	300	376
4	50	0.58	300	404

RESULT

The value of Stefan's constant, with the help of the experiment is

$$\sigma = 5.17 \times 10^{-5} \text{ ergs cm}^{-2} \text{ sec}^{-1} \text{ K}^{-4}$$

$$\text{Standard Value} = 5.67 \times 10^{-5} \text{ ergs cm}^{-2} \text{ sec}^{-1} \text{ K}^{-4}$$

$$\text{Percentage error} = 8.8\%$$

PRECAUTIONS

1. The readings of thermometers are taken only as stated above.
2. All the fittings should be airtight.

PARTS LIST

- | | | |
|----|---------------------|--------------------------|
| 1. | MAIN UNIT | <input type="checkbox"/> |
| 2. | HEATER PLATES | <input type="checkbox"/> |
| 3. | THERMOMETER (3 nos) | <input type="checkbox"/> |
| 4. | CONNECTING WIRES | <input type="checkbox"/> |
| 5. | INSTRUCTION MANUAL | <input type="checkbox"/> |