

Benedict's reagent

Benedict's reagent (often called **Benedict's qualitative solution** or **Benedict's solution**) is a chemical reagent and complex mixture of sodium carbonate, sodium citrate, and copper(II) sulfate pentahydrate.^[1] It is often used in place of Fehling's solution to detect the presence of reducing sugars. The presence of other reducing substances also gives a positive result.^[2] Such tests that use this reagent are called the **Benedict's tests**. A positive test with Benedict's reagent is shown by a color change from clear blue to brick-red with a precipitate.

Benedict's test



Positive Benedict's test

Classification	Colorimetric method
Analytes	Reducing sugars

Generally, Benedict's test detects the presence of aldehydes, alpha-hydroxy-ketones, and hemiacetals, including those that occur in certain ketoses. Thus, although the ketose fructose is not strictly a reducing sugar, it is an alpha-hydroxy-ketone and gives a positive test because the base in the reagent converts it into the aldoses glucose and mannose. Oxidation of the reducing sugar by the cupric (Cu^{2+}) complex of the reagent produces a cuprous (Cu^+), which precipitates as insoluble red copper(I) oxide (Cu_2O).^[3]

It is named after American chemist Stanley Rossiter Benedict.^[4]

Contents

Composition and preparation

Organic analysis

Quantitative reagent

Net reaction

See also

References

Composition and preparation

Benedict's reagent is a deep-blue aqueous solution. Each litre contains:^[4]

- 17.3g copper sulfate
- 173g sodium citrate
- 100g anhydrous sodium carbonate or, equivalently, 270g sodium carbonate decahydrate

Separate solutions of the reagents are made. The sodium carbonate and sodium citrate are mixed first, and then the copper sulfate is added slowly with constant stirring.

Sodium citrate acts as a complexing agent which keeps Cu^{2+} in solution, since it would otherwise precipitate. Sodium carbonate serves to keep the solution alkaline. In the presence of mild reducing agents, the copper(II) ion is reduced to copper(I), which precipitates in the alkaline conditions as very

conspicuous red copper(I) oxide.

Organic analysis

To test for the presence of monosaccharides and reducing disaccharide sugars in food, the food sample is dissolved in water and a small amount of Benedict's reagent is added. During a water bath, which is usually 4–10 minutes, the solution should progress through the colors of blue (with no reducing sugar present), orange, yellow, green, red, and then brick red precipitate or brown (if a high concentration of reducing sugar is present). A color change would signify the presence of a reducing sugar.^[2]

Experiment	Observation	Inference
Substance in water + 3 mL Benedict's solution, then boil for few minutes and allow to cool.	Red, green, or yellow precipitate is obtained	Reducing sugar, such as glucose, is present
Substance in water + 3 mL Benedict's solution, then boil for few minutes and allow to cool.	Solution remains clear or is a little blue	Reducing sugar is not present

The common disaccharides lactose and maltose are directly detected by Benedict's reagent because each contains a glucose with a free reducing aldehyde moiety after isomerization.

Sucrose (table sugar) contains two sugars (fructose and glucose) joined by their glycosidic bond in such a way as to prevent the glucose undergoing isomerization to an aldehyde, or fructose to alpha-hydroxy-ketone form. Sucrose is thus a non-reducing sugar which does not react with Benedict's reagent. However, sucrose indirectly produces a positive result with Benedict's reagent if heated with dilute hydrochloric acid prior to the test, although it is modified during this treatment as the acidic conditions and heat break the glycosidic bond in sucrose through hydrolysis. The products of sucrose decomposition are glucose and fructose, both of which can be detected by Benedict's reagent as described above.

Starches do not react or react very poorly with Benedict's reagent due to the relatively small number of reducing sugar moieties which occur only at the ends of carbohydrate chains. Other carbohydrates which produce a negative result include inositol.

Benedict's reagent can also be used to test for the presence of glucose in urine, elevated levels of which is known as glucosuria. Glucosuria can be indicative of diabetes mellitus, but Benedict's test is not recommended or used for diagnosis of the aforementioned condition. This is due to the possibility of a reaction in which the presence of other reducing substances such as ascorbic acid, drugs (levodopa, contrast used in radiological procedures) and homogentisic acid (alkaptonuria) creates a false positive.


As color of the obtained precipitate can be used to infer the quantity of sugar present in the solution, the test is semi-quantitative. A greenish precipitate indicates about 0.5 g% concentration; yellow precipitate indicates 1 g% concentration; orange indicates 1.5 g% concentration; and red indicates 2 g% or higher concentration.

Quantitative reagent

Benedict's quantitative reagent contains potassium thiocyanate and is used to quantitatively determine the concentration of reducing sugars.^[2] This solution forms a copper thiocyanate precipitate which is white and can be used in titration. The titration should be repeated with 1% glucose solution instead of the sample for calibration.

Net reaction

The net reaction between an aldehyde (or an alpha-hydroxy-ketone) and the copper(II) ions in Benedict's solution may be written as:



The hydroxide ions in the equation forms when sodium carbonate dissolves in water. With the citrate included, the reaction becomes:



See also

Other oxidizing reagents

- Dextrose equivalent
- Fehling's solution
- Tollens' reagent

Other reducing reagents

- Jones reductor
- Walden reductor

References

1. Robert D. Simoni; Robert L. Hill & Martha Vaughan (2002). "Benedict's Solution, a Reagent for Measuring Reducing Sugars: the Clinical Chemistry of Stanley R. Benedict" ([https://www.jbc.org/article/S0021-9258\(19\)61050-1/fulltext](https://www.jbc.org/article/S0021-9258(19)61050-1/fulltext)). *J. Biol. Chem.* **277** (16): 10–11. doi:10.1016/S0021-9258(19)61050-1 (<https://doi.org/10.1016%2FS0021-9258%2819%2961050-1>).
2. Collins Edexcel International GCSEBiology, Student Book (ISBN 978-0-00-745000-8) p.42-43
3. "Carbohydrates - Benedict's Test" (<http://dept.harpercollege.edu/chemistry/chm/100/dgodambe/thedisk/carbo/bened/benedict.htm>). *dept.harpercollege.edu*. Retrieved 2020-03-08.
4. Benedict, S. R. (1 January 1909). "A Reagent For the Detection of Reducing Sugars" ([https://www.jbc.org/article/S0021-9258\(18\)91645-5/pdf](https://www.jbc.org/article/S0021-9258(18)91645-5/pdf)). *J. Biol. Chem.* **5** (6): 485–487. doi:10.1016/S0021-9258(18)91645-5 (<https://doi.org/10.1016%2FS0021-9258%2818%2991645-5>).
5. <https://microbenotes.com/benedicts-test/>

Retrieved from "https://en.wikipedia.org/w/index.php?title=Benedict%27s_reagent&oldid=1026265177"

This page was last edited on 1 June 2021, at 08:33 (UTC).

Text is available under the Creative Commons Attribution-ShareAlike License; additional terms may apply. By using this site, you agree to the Terms of Use and Privacy Policy. Wikipedia® is a registered trademark of the Wikimedia Foundation, Inc., a non-profit organization.