

## Exp. 12. Coefficient of viscosity of water by Poiseuille's method

### Object

To determine the coefficient of viscosity of water by the Poiseuille's method.

### Apparatus

Large glass jar with a hole at the bottom, meter scale, capillary tube, rubber, stop-watch, spirit level and a beaker.

### Theory

A liquid opposes relative motion between its different layers due to its internal friction called viscosity, and, therefore, force has to be applied to maintain its steady motion. When a liquid flows through a narrow tube under ordinary pressure, its layer in contact with the tube is stationary and in going towards the axis of the tube, its velocity increases, and it is same at all points at the same distance from the axis of the tube.




Fig. 12.1

Let us consider the flow of liquid between two adjacent cylindrical layers coaxial with the tube and of radii  $r$  and  $r + dr$  (fig. 12.1). According to Newton's law,

Viscous force  $F = -\eta \times \text{area} \times \text{velocity gradient}$

$$\text{or, } F = -\eta \times 2\pi r l \times \frac{dv}{dr}$$

where  $\eta$  is the coefficient of viscosity and  $l$  is the length of the tube.

For steady motion,

$$F = p\pi r^2$$

where  $p$  is the pressure difference between the ends of the tube.

$$\therefore p\pi r^2 = -\eta \times 2\pi r l \times \frac{dv}{dr}$$

$$\text{or } dv = -(p/2\eta l) r dr$$

$$\text{Integrating, } v = -(p/2\eta l) \frac{r^2}{2} + A$$

[A is the constant of integration]

Now, near the wall,  $v = 0$  and  $r = a$  = the radius of the tube.

$$\therefore A = (pa^2)/(4\eta l)$$



[B] Readings for the determination of  $h$  and  $Q$

Height of the capillary tube above the experimental table,  $h_3 = \dots$  cm

S. No.	* $h_1$ cm	* $h_2$ cm	$h =$ cm	$V =$ Volume of water collected in c.c.	$t =$ Time taken in sec.	Vol. flowing per sec. $\left(Q = \frac{V}{t}\right)$	$\frac{h}{Q}$
1.							
2.							
3.							

\*While using constant head water tank, these two columns must be deleted.

[C] (i) Length of the capillary tube = ... cm

(ii) Room temperature = ... °C

### Calculations

Mean radius of the tube ( $a$ ) = ... cm

Now,

$$\eta = \frac{\pi(hpg)a^4}{8Ql}$$

$$= \frac{\pi ga^4}{8l} \left( \frac{h}{Q} \right)$$

[since  $p \approx 1$ ]

= ... dyne per cm<sup>2</sup> per unit velocity gradient (or poise)

% error :

### Result

The coefficient of viscosity of water ... °C (correct to significant figures) = ... poise.

### Precautions

- The capillary tube should be placed horizontally.
- The capillary tube should have a uniform bore about 0.5 mm in diameter, otherwise the flow becomes turbulent.
- Care should be taken to avoid back lash error in the microscope while taking the reading for diameter of the tube.
- Too much water should not be collected, otherwise  $h_1$  and  $h_2$  will differ much and the average pressure difference between the two ends of the capillary tube will not be constant.
- The pinch cock should be opened completely while determining the rate of water through the capillary tube.

### Important points about the experiment

There are two main sources of error in the above experiment [a] part of the thrust, due to the pressure difference between the two ends of the capillary tube, imparts kinetic energy to the liquid and the whole of which is not used in  $m$  overcoming the viscous resistance of the liquid. This may be corrected replacing  $p$

by  $p - \frac{Q^2 \rho}{\pi^2 a^4}$ , [b] the motion of the liquid, where it enters the capillary tube, is accelerated with the result that the velocity of flow is not uniform for the short length of the tube. This is eliminated by replacing  $l$  by  $[l + 1.64a]$ . Thus  $\eta$  becomes

$$\eta = \frac{\pi p a^4}{8Q(l + 1.64a)} = \frac{Qp}{8\pi(l + 1.64a)}$$

*viscosity is a quantity expressing the magnitude of internal friction in a fluid, as measured by the force per unit area resisting uniform flow.*

*The degree to which a fluid resists flow under an applied force, measured by tangential friction per unit area divided by the velocity gradient under the conditions of streamline flow.*

### QUESTIONS

- Define viscosity and coefficient of viscosity.
- What are the units and dimensions of coefficient of viscosity and how is it affected by temperature and pressure? *SI → poiseuille, CGS → poise, [ML<sup>-1</sup>T<sup>-1</sup>]. increases with increase in pressure as reported to temp.*
- Distinguish between streamline and turbulent motions.
- What are the conditions on which the validity of Poiseuille's formula depends and how can they be fulfilled in practice? What are its limitations?
- Why should you take a glass jar of large diameter and a capillary tube of 0.5 mm diameter? Is this method suited to determine the viscosity of all liquids?
- What limits the accuracy of the result in this experiment?
- Would you prefer to use a long tube or a small tube?
- Why is the capillary tube kept horizontal? *To keep the flow streamline.*

### Exp. 13A. Coefficient of viscosity of Liquid

#### Object

To determine the coefficient of viscosity of a liquid by rotating viscometer (Searle's apparatus).

#### Apparatus

Searle's apparatus, lubricating oil, vernier callipers, experimental liquid, a stop-watch and weight box.

The apparatus consists of two coaxial metal cylinders A and B, the space between which contains the experimental liquid (fig. 13.1). The inner solid cylinder A is fixed to an axle E which is pivoted freely at its ends. The axle can be attached to